Omegamatrix-Methode

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öffentlich

OM04-000a_v1.0 gültig ab 26-Mrz-2025

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öffentlich

OM04-000a_v1.0 gültig ab 26-Mrz-2025

Influencing and interfering factors / potential residual risk

Like any analytical method, our Omegamatrix Method is also subject to influencing and interfering factors that can affect the measurement result. Specific measures have therefore been taken to minimize these factors and ensure the accuracy and reliability of the results.

1. Mix-up of sample types (EDTA blood vs. dried blood):

The Omegamatrix Method can be performed with venous blood (EDTA blood) and dried capillary blood (dried blood). With EDTA blood, the fatty acids in the erythrocytes can be analyzed separately from the plasma, whereas dried blood contains the total blood matrix. This leads to systematic deviations, which are compensated for by a sliding correction factor. Mixing up the samples can nevertheless lead to erroneous results. Therefore, samples must be clearly labelled and processed accordingly.

2. Improper storage of EDTA blood samples:

The storage of EDTA blood samples below 0°C has an influence on the measurement result. Therefore, EDTA blood samples must not be stored below 0°C before and during return to the Omegamatrix laboratory. The reason for this is the freezing and thawing of the sample and the associated bursting of the erythrocytes. This is noticeable after centrifugation of the EDTA blood samples due to a lack of phase formation. Affected samples can therefore be easily identified and are not used for analysis. These EDTA blood samples are discarded and new sample material is requested from the sender.

3. Insufficient fasting before collection of capillary blood:

Note: This does not apply to the collection of EDTA blood.

Food intake prior to the collection of capillary blood influences the measurement result, as it is not possible to separate the erythrocytes from the rest of the whole blood in the laboratory when preparing dried blood samples (dried capillary blood). As the fatty acids in the erythrocytes are measured, increased fatty acids in the plasma can falsify the result. In order to minimize this influence and reduce the biological variability of the measurement, it is therefore essential that capillary blood is taken from the fingertip in a fasting state, i.e. after fasting for at least 10 hours.

Potential residual risk:

Despite the consideration of all known interfering factors and compliance with pre-analytical specifications, a certain residual risk always remains when determining parameters. A 100% certainty that unknown or unforeseeable interfering factors will not influence the result cannot be guaranteed. Therefore, this residual risk should always be included in a possible error analysis.